



# Medical Device Alert

Ref: MDA/2012/062 Issued: 11 September 2012 at 15:00

## Device

Silicone tracheostomy tubes

Various manufacturers

Various models

## Problem

The MHRA is aware of several incidents where the eyelets of tracheostomy tubes appear to have been damaged by the holder. This led to dislodgement of the tube with the subsequent risk of airway loss.

## Action by

All staff responsible for the care of patients with tracheostomy tubes.

## CAS deadlines

Action underway: 25 September 2012

Action complete: 09 October 2012

**Note: These deadlines are for systems to be in place to take actions.**

## Action

- Check the instructions for use for both the tracheostomy tube and the holder to ensure that they are compatible.
- Be aware that some manufacturers have recently updated their instructions for use to include a warning to users not to use their tracheostomy tubes with sharp-edged holders, for example some types of Velcro or metal edged holders.
- Ensure that the instructions for use are followed for both the tracheostomy tube and the holder.

## Device

Tracheostomy tubes are held in position using holders (also known as neck ties). The holders are attached by being threaded through the eyelets on the flanges of the tracheostomy tube.

## Distribution

This MDA has been sent to:

- NHS trusts in England (Chief Executives)
- Care Quality Commission (CQC) (Headquarters) for information
- HSC trusts in Northern Ireland (Chief Executives)
- NHS boards in Scotland (Equipment Co-ordinators)
- Local authorities in Scotland (Equipment Co-ordinators)
- NHS boards and trusts in Wales (Chief Executives)
- Primary care trusts in England (Chief Executives)

### Onward distribution

Please bring this notice to the attention of relevant employees in your establishment. Below is a suggested list of recipients.

#### Trusts

CAS and SABS (NI) liaison officers for onward distribution to all relevant staff including:

- Adult intensive care units
- All clinical departments
- All clinical staff
- All wards
- Ambulance services directors
- Ambulance staff
- Anaesthetists
- Clinical governance leads
- Day surgery units
- Health and safety managers
- Hospital at home units
- Medical directors
- Medical libraries
- Nursing executive directors
- Outpatient departments
- Paediatric intensive care units
- Purchasing managers
- Risk managers
- Supplies managers
- Theatres

#### Primary care trusts

CAS liaison officers for onward distribution to all relevant staff including:

- Community children's nurses
- Community hospitals
- Community nurses
- District nurses
- General practitioners
- Health visitors
- Minor injury units
- NHS walk-in centres
- Practice nurses
- School nurses
- Walk-in centres

#### Social services

Liaison officers for onward distribution to all relevant staff including:

- In-house residential care homes

### Independent distribution

#### Establishments registered with the Care Quality Commission (CQC) (England only)

This alert should be read by:

- Care homes providing nursing care
- Care homes providing personal care
- Hospices
- Hospitals in the independent sector
- Independent treatment centres
- Private medical practitioners

Issued: 11 September 2012 at 15:00

Ref: **MDA/2012/062**

### **Establishments registered with OFSTED**

This alert should be read by:

- Educational establishments with beds for children
- Residential special schools

Please note: CQC and OFSTED do not distribute these alerts. Independent healthcare providers and social care providers can sign up to receive MDAs directly from the Department of Health's Central Alerting System (CAS) by sending an email to: [safetyalerts@dh.gsi.gov.uk](mailto:safetyalerts@dh.gsi.gov.uk) and requesting this facility.

## **England**

If you are in England, please send enquiries about this notice to the MHRA, quoting reference number **MDA/2012/062** or **2012/001/018/401/003**

### **Technical aspects**

Louise Mulroy and Gica Leclerc  
Medicines & Healthcare products Regulatory Agency  
Floor 4  
151 Buckingham Palace Road  
London SW1W 9SZ  
Tel: 020 3080 7344 and 6610  
Fax: 020 8754 3965  
Email: [louise.mulroy@mhra.gsi.gov.uk](mailto:louise.mulroy@mhra.gsi.gov.uk)  
[gica.leclerc@mhra.gsi.gov.uk](mailto:gica.leclerc@mhra.gsi.gov.uk)

### **Clinical aspects**

Jonathan Plumb  
Medicines & Healthcare products Regulatory Agency  
Floor 4  
151 Buckingham Palace Road  
London SW1W 9SZ  
Tel: 020 3080 6648  
Fax: 020 8754 3965  
Email: [jonathan.plumb@mhra.gsi.gov.uk](mailto:jonathan.plumb@mhra.gsi.gov.uk)

## **How to report adverse incidents**

Please report via our website <http://www.mhra.gov.uk>

Further information about **CAS** can be found at <https://www.cas.dh.gov.uk/Home.aspx>

## **Northern Ireland**

Alerts in Northern Ireland will continue to be distributed via the NI SABS system.

Enquiries and adverse incident reports in Northern Ireland should be addressed to:

Northern Ireland Adverse Incident Centre  
Health Estates Investment Group  
Room 17  
Annex 6  
Castle Buildings  
Stormont Estate  
Dundonald BT4 3SQ  
Tel: 02890 523 704  
Fax: 02890 523 900  
Email: [NIAIC@dhsspsni.gov.uk](mailto:NIAIC@dhsspsni.gov.uk)  
<http://www.dhsspsni.gov.uk/index/hea/niaic.htm>

## **How to report adverse incidents in Northern Ireland**

Please report directly to NIAIC, further information can be found on our website <http://www.dhsspsni.gov.uk/niaic>  
Further information about **SABS** can be found at <http://sabs.dhsspsni.gov.uk/>

## **Scotland**

Enquiries and adverse incident reports in Scotland should be addressed to:

Incident Reporting and Investigation Centre

Health Facilities Scotland

NHS National Services Scotland

Gyle Square

1 South Gyle Crescent

Edinburgh EH12 9EB

Tel: 0131 275 7575 Fax: 0131 314 0722

Email: [nss.irc@nhs.net](mailto:nss.irc@nhs.net)

<http://www.hfs.scot.nhs.uk/online-services/incident-reporting-and-investigation-centre-irc/>

## **Wales**

Enquiries in Wales should be addressed to:

Improving Patient Safety Team

Medical Directorate

Welsh Government

Cathays Park

Cardiff CF10 3NQ

Tel: 029 2082 3922

Email: [Haz-Aic@wales.gsi.gov.uk](mailto:Haz-Aic@wales.gsi.gov.uk)